











Analytical columns with NUCLEOSIL® C₁₈ phases



Length →	100 mm	125 mm	150 mm	250 mm	Guard columns
NUCLEOSIL® 100-5 C₁₈ particle size 5 µm, pore size 100 Å, 15 % C					
EC columns					
	2 mm ID	720002.20		720014.20	721602.30
	3 mm ID	720002.30		720014.30	721602.30
	4 mm ID	720002.40		720014.40	721602.40
	4.6 mm ID	720002.46	720120.46	720014.46	721602.40
ChromCart® cartridges					
	2 mm ID	721622.20		721662.20	721602.30
	3 mm ID	721622.30		721662.30	721602.30
	4 mm ID	721622.40		721662.40	721602.40
	4.6 mm ID	721622.46	721642.46	721662.46	721602.40
Microbore columns					
	1 mm ID	717028.10	717019.10	717010.10	717002.10
NUCLEOSIL® 100-7 C₁₈ particle size 7 µm, pore size 100 Å, 15 % C					
EC columns					
	4 mm ID			720018.40	
	4.6 mm ID			720018.46	
ChromCart® cartridges					
	3 mm ID	721878.30		721609.30	
	4 mm ID	721878.40		721609.40	
	4.6 mm ID			721609.46	
NUCLEOSIL® 100-10 C₁₈ particle size 10 µm, pore size 100 Å, 15 % C					
EC columns					
	4 mm ID			720023.40	
	4.6 mm ID			720023.46	
ChromCart® cartridges					
	3 mm ID	721681.30		721689.30	
	4 mm ID	721681.40		721689.40	
	4.6 mm ID			721689.46	
NUCLEOSIL® 120-3 C₁₈ particle size 3 µm, pore size 120 Å, 11 % C					
EC columns					
	2 mm ID	720040.20		720055.20	721606.30
	3 mm ID	720040.30		720055.30	721606.30
	4 mm ID	720040.40		720055.40	721606.40
	4.6 mm ID	720040.46	720740.46	720055.46	721606.40
ChromCart® cartridges					
	2 mm ID	721626.20		721666.20	721606.30
	3 mm ID	721626.30		721666.30	721606.30
	4 mm ID	721626.40		721666.40	721606.40
	4.6 mm ID	721626.46	721646.46	721666.46	721606.40
Microbore columns					
	1 mm ID	717031.10	717022.10	717013.10	717005.10

8 mm ChromCart® guard column cartridges in packs of 3, all other columns in packs of 1.

As guard columns for EC columns use ChromCart® guard column cartridges with guard column adaptor EC (Cat. No. 721359). ChromCart® columns require the CC connecting kit (Cat. No. 721690).

On request, Microbore columns are also available in lengths of 40, 60, 200 and 300 mm and with 0.05, 0.075, 0.1, 0.15, 0.3, 0.4, 0.5, 0.75 and 1.5 mm ID. Guard columns for Microbore columns on request.